
Leaflet

LF-EXT-057
Version 10.1

CHROMagar™ mSuperCARBA™

Click below:

CHR  Magar™
The Chromogenic Media Pioneer

CHROMagar™ mSuperCARBA™

**For detection and isolation of Carbapenemase-Resistant
Enterobacteriaceae (CRE)**





Plate Reading

-  • CPE *E. coli*
→ dark pink to reddish
-  • CPE *Klebsiella, Enterobacter, Citrobacter*
→ metallic blue
-  • CPO* *Pseudomonas*
→ translucent cream to blue
-  • CPO* *Acinetobacter*
→ cream, opaque

*CPO : carbapenemase producer organism

For detection and isolation of Carbapenemase-resistant *Enterobacteriaceae* (CRE)

Background

CDC: «Carbapenem-resistant *Enterobacteriaceae* (CRE) are usually resistant to all β -lactam agents as well as most other classes of antimicrobial agents. The treatment options for patients infected with CRE are very limited. Healthcare-associated outbreaks of CRE have been reported. Patients colonized with CRE are thought to be a source of transmission in the healthcare setting. Identifying patients who are colonized with CRE and placing these patients in isolation precautions may be an important step in preventing transmission».

CHROMagar™ launched in 2007 the first chromogenic medium for the detection of carbapenem-resistant bacteria, particularly targeting KPC-enzymes. Since then, many other carbapenemases have been spreading around the world and therefore there was a need to address today the difficult detection of low level carbapenemases.

Alain Rambach and Patrice Nordmann have joined their efforts to develop a highly sensitive chromogenic medium, CHROMagar™ mSuperCARBA™, the new generation of chromogenic media that reaches unprecedented performances: detection of a large variety of carbapenemases KPC, NDM, VIM, IMP, OXA...with an impressive limit of detection (10 CFU/mL) even for weakly expressed carbapenemases like OXA-48, while maintaining a high level of selectivity.

Intended Use

CHROMagar™ mSuperCARBA™ is a selective and differential chromogenic culture medium, intended for use in the qualitative direct detection of gastrointestinal colonization with carbapenem-resistant *Enterobacteria* (CRE), including OXA-48 producers, to aid in the prevention and control of CRE in healthcare settings. The test is performed with rectal swab and stools from patients to screen for CRE colonization. Results can be interpreted after 18-24 h of aerobic incubation at 35-37 °C.

CHROMagar™ mSuperCARBA™ is not intended to diagnose CRE infection nor to guide nor monitor treatment for infections. A lack of growth or the absence of colonies on CHROMagar™ mSuperCARBA™ does not preclude the presence of CRE. Further identification, susceptibility testing, and epidemiological typing is needed on suspect colonies.

Medium Performance

- 1 HIGHLY SENSITIVE**
Most carbapenemases (including OXA-48 and OXA-48 like) detected after overnight incubation.
- 2 IMPRESSIVE LIMIT OF DETECTION**
10 CFU/mL
- 3 HIGHLY SELECTIVE AND SPECIFIC**
Not only beta-lactam susceptible bacteria are inhibited, but also most of the ESBL and AmpC hyperproducers are inhibited, rendering the tool very specific for carbapenemase-producing bacteria.

	Analytical data *	Clinical Data **
	CHROMagar™ mSuperCARBA™	
Sensitivity	100 %	100 %
Specificity	71 %	100 %

* Data obtained after 24 h incubation at 37 °C in aerobic conditions in the study «Amélioration de la détection des Entérobactéries Productrices de Carbapénémase (EPC)». Dos Santos *et al.* RICA 2017.

** Data obtained after 24 h incubation at 35 °C with 211 rectal swabs from the study «CHROMagar™ mSuperCARBA: performance in carbapenem-resistant *Enterobacteriaceae* isolates characterized at molecular level and routine surveillance rectal swabs specimens». García-Fernández *et al.* 2017. *Diagn. Microbiol. Infect. Dis.*

Medium Description

Powder Base	Total 42.5 g/L Agar 15.0 Peptones 20.0 Salt 5.0 Chromogenic and selective mix 0.8 Growth factors 1.7 Storage at 15/30 °C - pH: 7.2 +/- 0.2 Shelf Life > 18 months	
2 Supplements (included in the pack)	1 st : Liquid form.....2 mL/L Storage at 15/30 °C Shelf Life > 18 months	2 nd : Powder..... 0.25 g/L Storage at 2/8 °C Shelf Life > 18 months

Usual Samples	Rectal swabs and stools.
Procedure	Direct Streaking, Incubation 18-24 h at 35-37 °C Aerobic conditions.

Scientific Publications on this product: available on www.CHROMagar.com
Please read carefully the instructions for use (IFU document) available on www.CHROMagar.com

Order References

Please use these product references when contacting your local distributor:

5000 mL pack SC172 (included in this reference: powder base SC172) + suppl. SC172(S1) + SC172(S2))

25 L pack SC173-25

(included in this reference: powder base SC173-25) + suppl. SC173-25(S1) + SC173-25(S2))

Manufacturer: CHROMagar
4 place du 18 juin 1940 75006 Paris - France
Email: CHROMagar@CHROMagar.com
Website: www.CHROMagar.com

Find your nearest distributor on
www.CHROMagar.com/contact

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