

ahn myTip®

Reload System RL



Eco-friendly and cost-effective tip management solution – each reload tower accommodates multiple insert plates of refill tips. This minimises plastic waste **by approximately 70% per pack** and reuses the current AHN myTip® RB rack boxes by simply replacing empty plates with full ones from the AHN myTip® RL reload towers

Secure fit of pipette tips – pipette tips sit stable and secure in the insert plate, enabling easy handling with multi channel pipettes

Secure and comfortable handling – AHN myTip® RL insert plate securely snaps into the rackbox, but can be easily removed again. To facilitate the handling for users with small hands, the tabs of the AHN myTip® RL insert plate are located on the longer side

Easy identification of pipette tips – colour-coded racks allows quick and easy identification of the pipette tip, and therefore a suitable pipette

Space-saving – small footprint of the AHN myTip® RL reload towers means there is more space both on the benchtop and in the storage room



metal free



1

Open the refill by sliding it towards the top

2

Prepare an empty rack. Remove its original grid

3

Pick up the first grid of refill tips and insert into the rack. The racks can be reused

4

Once the tips are properly installed in the rack, you may close the hinged lid

ahn



+49(0)3631/65242-0

+49(0)3631/65242-90

Uthleber Weg 14 | 99734 Nordhausen | Germany



info@cappahn.com






www.ahn-bio.de



Ordering information

Reload System (RL)

Volume	Length	Colour Tip	Colour Code	Sterilised	Pack Type	Sales Unit	Cat. No.
 AHN myTip® RL 10 µL	31.20mm	Clear	■		reload	case / 10 x 10 x 96 tips	1-002-C1-R
 AHN myTip® RL 10/20 µL XL	45.80mm	Clear	■		reload	case / 10 x 10 x 96 tips	1-062-C1-R
 AHN myTip® RL 200 µL	53.40mm	Clear	■		reload	case / 10 x 10 x 96 tips	1-121-C1-R
 AHN myTip® RL 300 µL	59.35mm	Clear	■		reload	case / 10 x 8 x 96 tips	1-127-C2-R
 AHN myTip® RL 1000 µL	84.30mm	Clear	■		reload	case / 10 x 5 x 96 tips	1-202-C3-R

ahn myTip®

Rack Box System RB

(without Tips)

Effortless and secure handling - easy one-handed opening and closing of rack box / secure rack box closing / rack boxes can be stacked for easy storage

Universal fit - suitable for multi channel pipettes (standard 9mm SBS footprint)

Easy identification of pipette tip - colour-coded insert plates allow for quick identification of the pipette tip volume






Rack boxes available in small & large size / racks in 8 x 12 format



metal free

Ordering information

Rack Box System without Tips (RB)

Volume	Colour	Sales Unit	Cat. No.
AHN myTip® RB for tips 10 µL/10 µL XL 	clear black (RB) / light grey (I)*	case / 5 x 10 rack boxes	4-500-C5-E
AHN myTip® RB for tips 200 µL/300 µL 	clear black (RB) / yellow (I)*	case / 5 x 10 rack boxes	4-501-C5-E
AHN myTip® RB for tips 200 µL Gelloader 	clear black (RB) / yellow (I)*	case / 4 x 8 rack boxes	4-502-C4-E
AHN myTip® RB for tips 1000 µL 	clear black (RB) / blue (I)*	case / 4 x 8 rack boxes	4-503-C4-E
AHN myTip® RB for tips 1250 µL 	clear black (RB) / blue (I)*	case / 4 x 8 rack boxes	4-504-C4-E

*RB - Rack Box | I - Insert Plate

ahn



+49(0)3631/65242-0

+49(0)3631/65242-90

Uthleber Weg 14 | 99734 Nordhausen | Germany



info@cappahn.com



www.ahn-bio.de

