


# EU Thin Wall (q)PCR 24 well semi-skirted plate


## For Diagnostic & Clinical procedures

e x c e l l e n t p e r f o r m a n c e a n d h a n d l i n g

- Perfect fit to nearly ALL PCR & qPCR cyclers  
(see table next page)
- Fits ABI 7300 and 7500 Real Time Cyclers
- Designed to perform qPCR and PCR reactions
- Cutable
- Low evaporation due to clever design
- Easy labeling and orientation, can be color coded by “clicking” the grid onto the plate(s)
- Extreme Uniform wells lead to reproducible results
- Can be closed using any of the BIOplastics cap strips
- Can be placed, transported & stored in Tube Support grid, Station Support Insert or Work Rack S-96



light frosted surface enables signal enhancement in qPCR



orientation & “click” hole

Cut positions

### EU Semi Skirted Thin-wall 24 x 0.2 ml plate, 3 x 8

Can be used in most 0.2 ml block regular and Real-Time thermal cyclers  
(see table next page )

**Orderno Description**

B50240	24 x 0.2ml Plate, Natural, Box of 100 plates, Semi-skirted Thin-wall
B57801	EU Wide Optical 8-Cap Strip, Wide Indented Flat Cap
B75701	EU Optical 8-Cap strip
B69351	EU Pre-Post Tube Support Grid Wide, Natural, Re-useable

**Package Size**

Box, 100 plates
Natural, Bag of 120 Strips
Natural, Bag of 120 Strips
Natural, 8 Grids

**Price**