

Premium Benchtop Shaking Incubators





ZWYR-240

ZWYR-200D

Description:

This new premium benchtop shaking incubators provided by LABWIT offer conditions of time, broad temperature and orbital or reciprocal agitation capabilities for meeting comprehensive requirement of the growth of biological organisms. This P.I.D microprocessor controlled shakers have a number of advanced features that are really unique in bench top models.


This unique benchtop model with shaking diameter stepless adjustable from 1-50mm provides unmatched flexibility for meeting comprehensive application requirements. The sophisticated PID controller provides these units with choices of not only constant controlling one fixed temperature and speed, but also programmed controlling with a series of "ramps and soaks" segments on time basis.

For data recording, the units are equipped with a built-in thermo printer as standard, so that operators could check how the performance parameters vary according to specified process without having to be on site.

Features:

- 4.3" LCD 480x272 touch screen panel clearly indicates all parameters in one page display and makes it easy
 to change the operating parameters settings, even those multi set points under programmable mode
 intuitively.
- P.I.D microprocessor controller provides unmatched versatility by enable users to create personalized program (with up to 9 segments, with cycling) to automate changes to function parameters.
- Inner chamber & tray from high graded #304 electro polished stainless steel.
- Robust single eccentric counterbalanced drive in cast iron housing provides vibration and trouble free operation with high shaking speed up to 600rpm.
- The versatility of the driving system allows stroke continuous adjustment from 1 to 50 mm for optimum
 oxygenation and nutrient mixing and agitation in from small to large flasks.
- Quality brushless AC inductive motor supports reliable, quiet and smooth shaking movements up to top speed.
- Auto-defrosting system, ensuring long-terms non-frosting operation with low temperature setpoints.
- Non-volatile memory saves settings during a power outage and automatically restarts the unit after power is restored.
- Audible and visual alarms on over temperature and motor overheating, and automatically cut off the power
- Soft start and smooth acceleration when lid is closed after short opens.
- For data recording, the units are equipped with a built-in thermo printer as standard, so that operators
 could check how the performance parameters vary according to specified process without having to be on
 site.
- Timer up to 9999 mins under fixed value control mode.
- Predrilled standard tray included, without clamps.

Flask (ml)	50ml	100ml	250ml	500ml	750ml	1000ml	2000ml	P8010	P8012
ZWYR-200D	23	23	12	9	7	5	4	3	8
ZWYR-240	23	23	12	9	7	5	4	3	8
ZWFR-200	23	23	12	9	7	5	4	3	8





Specifications:

Model	ZWYR-200D ZWYR-240		ZWFR-200				
Control	P.I.D Microprocessor						
Control Mode	Fix Value or Program (up to 9 Segments)						
Control Panel	4.3" LCD Touch Screen						
Volume (L)	69						
Shaking Speed (rpm)	30-600 30-400		30-240				
Stroke (mm)	1-50 Steples	1-40 Stepless Adjustment					
Shaking Mode	Orb	Reciprocal					
Tray (mm) (WxD)	400	370	435x420				
Ambient Temperature(ºC)	5-25						
Temperature Range (ºC)	4-60						
Temperature Accuracy (ºC)	0.1						
Temperature Uniformity (ºC)	±1ºC@37°C						
Timer	1 to 9999 mins						
Inner Dimensions (mm) (WxDxH)	490x470x320	490x470x335	490x470x320				
Exterior Dimensions (mm) (WxDxH)	774x740x610						
Packing Dimensions (mm) (WxDxH)	844x810x770						
Net/Gross Weight (kg)	100/132						
Power (W)	700						
Electricity	220-240V 50/60 Hz						
Approval		CE, ISO	CE, ISO				
	Over-temperature Protection, Compressor Overload Protection, Electrical						
Security	Leakage Protection						
Additional	Built-in Printer, RS485 COM Kit (Option)						



Touch Screen Panel



Forced Air Convection



1-50mm Shaking Diameter

Adjustment

Part NO.	Description		Description	
ZWYR-240	ZWYR-240, 69L, Premium Orbital Shaking Incubator 4-60, 30-400rpm		T Clamp, S/S, for 150ml Flask, with Spring Retainer	
ZWYR-200D	ZWYR-200D,69L, Premium Orbital Shaking Incubator 4-60, 30-600rpm		T Clamp, S/S, for 300ml Flask, with Spring Retainer	
ZWFR-200	ZWFR-200,69L, Premium Reciprocal Shaking Incubator 4-60, 30-240rpm		T Clamp, S/S, for 500ml Flask, with Spring Retainer	
P5016	RS-485 Interface	P8005	T Clamp, S/S, for 750ml Flask, with Spring Retainer	
P6004	Tray for ZWYR-240/200D, ZWY-240/200D, S/S*	P8006	T Clamp, S/S, for 1000ml Flask, with Spring Retainer	
P6005	Tray for ZWFR-200, ZWF-200, S/S	P8007	T Clamp, S/S, for 2000ml Flask, with Spring Retainer	
P7003	Universal Tray for ZWYR-240/200D, ZWY-240/200D	P8010	Tube Rack S/S	
P7005	Universal Tray for ZWFR-200, ZWF-200	P8011	Tube Rack ABS Plastic	
P8001	T Clamp, S/S, for 50ml Flask, with Spring Retainer	P8012	Clamp for 96 Well Plate, S/S	
P8002	T Clamp, S/S, for 100ml Flask, with Spring Retainer	P8017	Sticky Mat, 20 x 20cm	
P8003	T Clamp, S/S, for 250ml Flask, with Spring Retainer		*S/S: Stainless Steel	







T Clamps P8017 P8010